



# RO-TAP Sieve Shaker

The RO-TAP Sieve Shaker is designed to provide Quality Control personnel with accurate and consistent particle analysis testing and repeatable results. The shaker has two operations: a horizontal, circular motion and a vertical, tapping motion. These operations allow material particles to stratify and "seek" critical openings in the test media. The Ro-TAP is designed for consistent, repeatable particle size testing.



## Models:

- RX-29 model = designed for testing with 8" diameter test sieves.
- RX-30 model = designed for testing with 12" diameter test sieves.

- Requires minimal maintenance
- Digital 99 minute timer/clock (accurate to 0.1 second)
- Heavy duty stationary unit
- Before shipment each unit is rest run and serialized

